



— Introducing —



**disposable
medical devices**

*better **health** and
a better **way of life***

www.umaflowindia.com

www.umaflowindia.com

*better **health** and
a better **way of life***



A revolution has taken place in the field of Medical devices with consistent innovation in technology.

Founded in 1993, **ROYAL SURGICARE PVT LTD.** Is a progressive and innovative company that deals with manufacturing & exporting medical devices offering the best quality services worldwide. Royal surgicare pvt ltd. boasts of an impressive State-of-the-art facility located in Santej, Gandhinagar, Gujarat INDIA.



INTRODUCTION

ROYAL SURGICARE PVT. LTD. offers a diversified full service manufacturing hub, with latest technology delivering superior quality medical devices.

We manufacture & export an extensive range of any of medical devices for infusion therapy, transfusion therapy, gastroenterology, urology, anesthesia, cardiology and surgery at our high tech facility situated in Santej, Gandhinagar, Gujarat INDIA.

In its drive to become the leading manufacturer & exporter in the field of Medical disposable, the company has developed quality and product System as per ISO 13485:2016 with all the products CE marked.


We deliver solutions with our expertise and unparalleled precision.

With royal surgicare pvt ltd.'s global footprint we are strategically positioned to help you realize your global initiatives.


Currently our global presence is felt in 35 countries across the world.

Come join the leaders in innovation and realize the dream of prioritizing Healthcare above all, "let's together make the world happy place to live".




T-7101
SAFE 'N' FLOW



- Chamber with air vent.
- Latex free "Y" injection site & luer lock.
- Available with or without needle.

T-7102
TOP


- Chamber without air vent.
- Tube latex injection site.
- Needle 21G.


T-7103
SEF-T


- 'T' shaped chamber without air vent.
- Bulb latex injection site.
- Needle 21G.

T-7104
1UP VENTED


- Chamber with air vent.
- Bulb latex / tube latex injection site.
- Needle 21G.
- Luer lock facility available on request.

T-7122
IV FLO



- Big chamber without air vent.
- Bulb latex injection site.
- Needle 21G.

T-7124
IV SET WITH 0.2 μ FILTER




- I.V. Set with 0.2 μ in-line filter.
- Chamber with air vent.
- Latex free "Y" injection site & luer lock.
- Also available with 60 drops / ml.




COMMON FEATURES

- Soft & kink free medical grade PVC tubing.
- Sharp and extra strong spike for ideal penetration in fluid container.
- Disk type fluid filter at the base of the chamber.
- Smooth roller clamp for controlling fluid rates and maintaining unrestricted flow.
- Standard length : 150 cm.
- Drop rate : 20 drops / ml.
- Customize packing available.
(Poly Bag, Paper Pouch, Soft Blister, Ribbon)

T-7105
"Y" MICRO



- Chamber with air vent.
- Latex free "Y" injection site & luer lock.
- Available with or without needle.

T-7106
MICRO DRIP


- Chamber without air vent.
- Bulb latex / tube latex injection site.
- Needle 23G.
- Also available with air vent chamber.

T-7107
BURETTE SET - 150 ML


- Measured volume set with 150ml burette capacity.
- Also available with latex free "Y" injection site & luer lock.
- Also available with air vent spike.

T-7108
BURETTE SET-110ML


- Measured volume set with 110ml burette capacity.
- Also available with latex free "Y" injection site & luer lock.
- Also available with air vent spike.

T-7113
INFUSION WITH DIAL FLOW



- Chamber with air vent.
- Dial flow regulator with flow rate control of IV fluid with range of 5 to 250 ml/hr.
- Latex free "Y" injection site & luer lock.
- Also available with 20 drops / ml.

COMMON FEATURES

- Soft & kink free medical grade PVC tubing
- Sharp and extra strong spike for ideal penetration in fluid container
- Disk type fluid filter at the base of the chamber
- Smooth flow regulator for controlling fluid rates and maintaining unrestricted flow.
- Standard length : 150 cm
- Drop rate : 60 drops/ml
- Customize packing available (Poly Bag, Paper Pouch, Soft Blister, Ribbon)



T-7109
BLOOD ADMINISTRATION SET - ECONOMY



- Single moulded chamber.
- Tube latex injection site.
- Also available with air vent on request.
- Needle 18G.

T-7119
BLOOD ADMINISTRATION SET "Y" INJECTION SITE



- Double moulded chamber.
- Latex free "Y" injection site & luer lock.
- Available with or without needle.

T-7110
BLOOD ADMINISTRATION SET - STANDARD



- Double moulded chamber.
- Bulb latex injection site.
- Needle 18G.

COMMON FEATURES

- Soft& kink free medical grade PVC tubing.
- Sharp and extra strong spike.
- Provided with in-built blood filter.
- Standard length : 150 cm.
- Smooth roller clamp for controlling and maintains unrestricted flow.
- Customize packing available.
(Poly Bag, Paper Pouch, Soft Blister, Ribbon)

T-7125
THREE WAY STOP COCK



- 360° rotation.
- Rotating male luer lock and two female luer lock connectors.
- For pressure up to 65psi.

T-7112
SCALP VEIN SET



- Designed to provide rapid venous access and greater patient comfort during infusion.
- Thin-walled beveled siliconized needle.
- Soft butterfly wings for easy handling.

Size	20G	21G	22G	23G	24G	25G
Needle Length	3/4"					

T-7114
HIGH PRESSURE EXTENSION LINE


- For pressure up to 800psi.
- With male & female connectors. • Kink free tubing.
- Available sizes : 10cm, 25cm, 50cm, 100cm, 150cm, 200cm.

T-7116/A
LOW PRESSURE EXTENSION LINE


- For pressure up to 65psi.
- With male & female connectors. • Kink free tubing.
- Available sizes : 10cm, 25cm, 50cm, 100cm, 150cm, 200cm.

T-7116
EXTENSION LINE WITH THREE WAY STOP COCK


- Three way stop cock at one end & male luer lock connector on other hand.
- Kink free tubing.
- Available sizes : 10cm, 25cm, 50cm, 100cm, 150cm, 200cm.

T-7118
IV CANNULA WITH INJECTION PORT & WINGS


- Beveled shaped needle for smooth penetration into the vein.
- Thin walled teflon catheter.
- Transparent flash back catheter for easy blood detection.

Size	14	16	18	20	22	24	26
Color	Orange	Grey	Green	Pink	Blue	Yellow	Violet
Length	45mm	45mm	45mm	32mm	25mm	19mm	19mm


T-7123
DISPOSABLE SYRINGES


- Sharp & clear graduations marked.
- Gaskets are made up of medical grade compound.
- Available with or without needle. • Slip luer tip.
- Size range : 2ml, 3ml, 5ml, 10ml, 20ml, 50ml.

T-7115
IV FLOW REGULATOR WITH EXTENSION SET


- Used to regulate the flow of iv fluid.
- Precise flow rate control of iv fluid with range of 5 to 250 ml/hr.
- Male luer lock at one end & female luer lock on the other hand.
- Latex free "Y" injection site.
- Standard length : 60cm.

T-7111
PERITONEAL DIALYSIS TRANSFUSION SET



- For administration of dialysis fluid.
- "Y" shaped set.

D-7701
A. V. FISTULA NEEDLE



- Tube length : 30cm.
- Sizes: 15G, 16G, 17G.
- Also available in twin pack.

D-7702
PERITONEAL DIALYSIS CATHETER KIT

- Stainless steel trocar is provided with triple angle sharpness to facilitate smooth penetration.
- Spare scalpel blade is supplied to provide the incision.
- The trocar and the tip of the catheter is perfectly matched to facilitate trauma free entry.
- Latex free flash ball is provided.
- Available sizes : adult and child.



GASTROENTEROLOGY


G-7308
UMBILICAL CATHETER

- Umbilical artery catheterization.
- Provides direct access to the arterial blood supply & allows accurate measurement of arterial blood pressure.
- Medical grade PVC tube with radio-opaque line.
- Color coded funnel connector length : 40 cm.
- Sizes : 3FG to 10FG.




G-7301
RYLE'S TUBE - ECONOMY



- Tube with radio - opaque line.
- Proximal end fitted with stopper.
- Cone shaped distal end with four lateral eyes.
- Length 105cm.
- Sizes : 6FG to 22FG.

G-7304
RYLE'S TUBE - UMALINE


- Tube with radio - opaque line.
- Proximal end fitted with closure.
- Cone shaped distal end with four lateral eyes.
- Length 105cm.
- Sizes : 6FG to 22FG.

G-7306
LEVIN'S TUBE


- Also known as gastro - duodenal feeding tube
- For gastro-intestinal feeding & aspiration.
- Medical grade PVC tube with radio-opaque line.
- Distilled end closed with four lateral eyes.
- Length : 125cm.
- Sizes : 6FG to 24FG.

G-7307
STOMACH TUBE


- For nasogastric introduction of nutrition.
- Medical grade PVC tube with radio - opaque line.
- Distilled end closed with two lateral eyes.
- Length : 80cm.
- Sizes : 6FG, to 24FG

G-7302
INFANT FEEDING TUBE

- Tube with radio opaque line.
- Distilled end is closed with two lateral eyes.
- Length - 53cm.
- Sizes : 4FG, to 10FG.



U-7202
UROFILL



- Rod type PVC thread hanger.
- Tube length 80cm.
- Top outlet with closure.

U-7203
ROYAL-29



- Milky type medical grade PVC urine bag cloth.
- Moulded hanger.
- Tube length - 100cm.
- Top outlet with closure.

U-7204
ROYAL-19



- Transparent medical grade PVC urine bag cloth.
- Moulded hanger with thread.
- Tube length - 90cm.
- Top outlet with closure.

U - 7212
ROYAL - 5

PUSH-Pull TYPE



- Transparent medical grade PVC urine bag cloth.
- Push-pull / T-type bottom outlet for emptying the urine bag.
- Tube length - 100cm.

U-7214
ROYAL - 9

Push-Pull Type



T-type

- Transparent medical grade PVC urine bag cloth.
- Tube length - 90cm.
- Also available with push-pull / T-type bottom outlet.

COMMON FEATURES

- Capacity - approx. 2000 ml.
- Soft & kink free medical grade PVC tubing all urine collecting bags provided with NVR. (non-return valve)
- Bags graduated in ml.
- Pvc tubes provided with inlet connector with cap.
- Customized urine bag capacity available on request.
- Customized tube length available on request.

U-7210
UROMEASURE URINE COLLECTING BAG WITH
MEASURED VOLUME METER



- Urine bag graduated in ml with 250 ml capacity of measured volume meter for accurate measurement of urine output.
- Provided with moulded hanger with bottom outlet. (push - pull / t-type).
- Tube length - 100 cm.
- Urine bag capacity approx. - 2000ml.
- PVC tube provided with inlet connector with cap.

U-7211
ROYAL CARE
LEG BAG SET



- Capacity approx. - 800ml.
- Provided with non-returnable valve.
- Available with leg straps.
- Tube length - 50cm.
- T-type bottom outlet drainage system.
- Bag graduated with ml.
- PVC tube provided with inlet connector with cap.

U-7208
URINE COLLECTING BAG - TINY



- For short-term urine collection.
- Suitable for both male & female.
- Capacity - 100ml.

U-7207
TUR-SET




- "Y" shape set for endoscopic irrigation during trans-urethral resection of prostate gland.
- Proximal end fitted with flexible latex tubing.

A-7408
URETHERAL CATHETER

- For short-term bladder catheterization through urethra.
- Length : 40cm.
- Sizes : 8FG to 20FG, K-90 & K-91.



U-7205
NELATON CATHETER



- For short-term bladder catheterization through urethra in males.
- Length - 40cm.
- Sizes : 6FG to 22FG.

U-7216
FEMALE CATHETER



- For short-term bladder catheterization through urethra in females.
- Length - 18cm.
- Sizes : 6FG to 22FG.

U-7215
RECTAL CATHETER



- For introduction of enema solution into rectum to release/aspirate rectal fluid.
- Length: 40cm.
- Sizes : 18FG to 32FG.

U-7206
MALE EXTERNAL CATHETER - URICATH



- External catheter used for urine incontinence for day & night use for male patients.
- Provided with self-adhesive.
- Sizes : small, medium, large, ex-large.

U-7213
FOLEY BALLOON CATHETER

- Made of silicon elastomer bonded with latex rubber.
- Distilled end is coned with super smooth eyes proportional to each other to provide non-traumatic intubation.
- Highly elastic balloon.
- Hard valve ensures easy inflation & deflation of balloon.
- Also available in 100% silicon.
- Sizes available
 - 2-way : 6FR, 8FR, 10FR, 12FR, 14FR, 16FR, 18FR, 20FR, 22FR, 24FR.
 - 3-way : 16FR, 18FR, 20FR, 22FR.




A-7406
GUEDEL AIRWAY



- Oro - pharyngeal airway.
- Rounded atraumatic edges.
- Sizes : 000, 00, 0, 1, 2, 3, 4, 5.

A-7404
TWIN BORE NASAL OXYGENSET



- Soft twin nasal prong tip.
- 200 cm star lumen pvc tube.
- Sizes : adult, child, infant, neonatal.

A-7403
NASAL OXYGEN CATHETER



- Open distilled end with multiple lateral eyes.
- Connector at proximal end.
- Length : 40cm. • Sizes : 8FG to 20FG.

A-7401
SUCTION CATHETER
STANDARD CONNECTOR



- To remove secretion from mouth, oropharynx, trachea & bronchial tube.
- Funnel type standard connector.
- Length : 53cm.
- Sizes : 5FG to 22FG.

A-7409
SUCTION CATHETER
THUMB CONTROL CONNECTOR



- To remove secretion from mouth, oropharynx, trachea & bronchial tube.
- Thumb control connector.
- Length : 53cm.
- Sizes : 5FG to 22FG.

**A-7405
OXYGEN FACEMASK**



- Clear and soft mask with adjustable elastic strip and integrated nose clip.
- 200cm long star - lumen medical grade PVC tubing. with connectors at both the end.
- Sizes : adult & paediatric.

**A-7413
HIGH CONCENTRATION - MASK**



- Specially designed for high concentration therapy.
- Clear and soft mask with 600ml reservoir bag.
- 200cm star - lumen medical grade pvc tubing.
- Sizes : adult & paediatric.

**A-7412
RESPIRATORY EXCISICER**



- Helps to develop respiratory fitness.
- Designed with three chambers to provide smooth exercises to patient.

**A-7415
NEBULIZER KIT**



- Kit consist of nebulizer mask with adjustable elastic strip & integrated nose clip, nebulizer chamber with nebulizer rate of approx. 3cc in 10 minutes and 200cm star-lumen medical grade PVC tubing.
- Sizes : adult & paediatric.

**A-7414
MULTI FLOW VENTURI MASK**




- Clear and soft mask with adjustable elastic strip and integrated nose clip.
- Provided with seven different oxygen concentration dilutes.
- 200cm star-lumen medical grade PVC tubing.
- Sizes : adult & paediatric.

**A-7423
CATHETER MOUNT**



- It connects patient ET to breathing circuit.
- Variants : standard & expandable.

**A-7411
ENDOTRACHEAL TUBE**

PLAIN

- Nasal/oral tip suitable for nasal / oral intubation.
- Sizes : 2.0 to 7.0mm.


CUFFED

- Provided with soft cuff at distal end.
- Non return valve with pilot cuff.
- Murphy eye at distal end.
- Sizes : 3.0 to 10.0mm.

COMMON FEATURES

- Radio-opaque line throughout the length to facilitate exact location of tube.
- Proximal end fitted with standard 15mm connector.
- Manufactured from thermo sensitive material readily conforming body curvature stable at body temperature.

**A-7416
BREATHING FILTER**
BVF FILTER

HMEF FILTER

HME FILTER


**A-7417
PLAIN**

**A-7418
SINGLE WATER TRAP**

**A-7419
DOUBLE WATER TRAP**

COMMON FEATURES

- It is the gas path way connecting patient's lung to the mechanical ventilator.
- 22mm corrugated, flexible and light weight tubing made of high-quality plastic.
- Sizes available: adult, pediatric.

A-7424
T-RECOVERY KIT


- "T" pc designed to use in recovery for the delivery of oxygen.
- Connects to 15mm et tube connector/ laryngeal mask.
- 200cm star-lumen medical grade PVC tubing.

A-7425
NEBULIZER "T" PC KIT


- Specially designed for nebulizing the patient orally.
- Kit consists of "T" shaped connector, mouth piece.
- Nebulizer chamber with rate of approx. 3cc, and 200cm star lumen medical grade PVC tubing.

A-7427
BAIN CIRCUIT

- The process includes patients inspiring fresh gas from the outer reservoir tube.
- 2lit. Antistatic re-breathing bag.
- 22M/15F – 15M fixed elbow and 3mtr. extra length oxygen tubing.
- Sizes : adult and pediatric.


A-7428
AMBU BAG (RESUSCITATOR BAG)

- The manual resuscitator is designed to provide effective breathing assistance.
- It consists of a silicone manual resuscitator, star lumen crush resistance oxygen tubing, reservoir bag, pressure limiting valve, and silicone mask.
- Sizes : adult, pediatric, infant.



**A-7429
SPINAL NEEDLE**


- Spinal anesthesia needle used during spinal anesthesia.
- Quincke point for smoother penetration.
- Needle length : 89mm.
- Sizes : 18G, 20G, 21G, 22G, 23G, 24G, 25G, 26G, 27G.

**A-7430
EPIDURAL NEEDLE**


- Touchy type epidural needle provided with clear depth marking for accurate insertion depth reading.
- Polished inner bevel edge needle minimizes the risk of Catheter shearing.
- Standard length : 80mm.
- Sizes : 16G, 17G, 18G, 19G, 20G.

**A-7431
EPIDURAL CATHETER**

- Epidural Catheter is specially designed for short term and long term continuous anesthesia and pain therapy.
- Luer Lock twist connector is supplied for easy and safe connection.
- Sizes : 16G & 18G.


**A-7432
EPIDURAL SET**

- Set consists of Epidural Needle, Epidural Catheter, LOR 10ml Syringe, Filter 0.22 μ with luer lock and catheter adaptor.



S-7501
CHEST DRAINAGE CATHETER
(INTRACOASTAL DRAINAGE CATHETER)



- Suitable for effective drainage after cardio thoracic surgery.
- DEHP free tube with radio opaque line with markings.
- Atraumatic and rounded open distal end.
- Proximal end is fitted with tapered connector.
- Sizes : 8, 10, 12, 14, 16, 18, 20, 22, 24, 28, 32, 36, 40.


CLOSED WOUND SUCTION UNIT

S-7502
UMAVAC SET



- Suitable for closed wound drainage under negative pressure post-operatively with the option to operate one or two catheters simultaneously.
- Graduated bellow with connector.
- Connecting tube with clamp and "Y" connector.
- Curved trocar with catheter.
- Bellow container capacity : 800ml and 400ml.
- Sizes : 10, 12, 14, 16, & 18FG.

S-7503
MINOR VAC SET




- Suitable for closed wound drainage under negative pressure.
- Used for minor surgery or plastic surgery.
- Bellow container capacity: 50ml.
- Provided with curved trocar with catheter.
- Sizes : 6 & 8FG.

S-7504 UMADRAIN
UNDER WATER SEAL DRAINAGE BAG



- It is used for collection of drainage fluid from thoracic cavity.
- Special hanger facility provided for easy carrying or hanging of the bag.
- PVC drainage bag with 1000ml capacity.
- Soft & thick medical grade PVC tubing.

S-7505
VACCUM SUCK SUCTION SET
YANKAUR SUCTION SET



- Suitable for pre-operative removal of secretion of body fluid.
- Yankaur handle is open with four lateral eyes at distal end.
- 2.5mtr. long tubing provided.
- Proximal end is fitted with funnel connector which can be fitted with any type of suction apparatus.
- Available with standard/crown tip.

S-7511
YANKAUR SUCTION HANDLE

CROWN TIP

STANDARD TIP

- Yankaur Handle is open with four lateral eyes.



S-7509 UMASEAL
UNDER WATER SEAL DRAINAGE SYSTEM

UMASEAL ADULT
2000ML
(DOUBLE CHAMBER)



UMASEAL MIDI
1200ML



UMASEAL CHILD
500ML



- Suitable for pleural drainage in conjunction with chest drainage catheters in cardio thoracic procedure.
- Graduation marked on bottle in ml to determine the amount of collection.
- Soft and kink free tubing ensures unrestricted flow.
- Stainless steel hanger provided for easy attachment with bed.
- Suction port is provided to connect with suction unit.

**G-7303
UMA ADK**


- Extra soft DEHP free drainage catheter with 2000ml collecting bag.
- Sizes : 16, 18, 20, 22, 24, 26, 28, 30, 32, 36.

**S-7510
REDON DRAIN CATHETER**


- Provided with radio opaque line with multiple smooth eyes.
- Marking on the tube at 19cm from distal end.
- Length: 50cm, 70cm, 120cm, 170cm.
- Sizes : 8, 10, 12, 14, 16 & 18fg

MISCELLANEOUS
**M-7605
CORRUGATED DRAINAGE SHEET**


- Multichannel drainage system for efficient drainage during operation.
- Provided with radio opaque line.
- Size : 25mm x 250mm.
- Customized size available.

**M-7601
INFANT MUCUS EXTRACTOR**


- Specially designed for aspiration of secretion from oropharynx in newly born babies.
- Chamber of mucus extractor is transparent, graduated with capacity of 25ml.

**M-7602
UMBILICAL CORD CLAMP**

- Suitable for clamping the umbilical cord after the child birth.
- Grooved clamping area to prepare the slipping of the umbilical cord clamp.



S-7505
KARMAN CANNULA(MTP CANNULA)



- Specially designed for aseptic medical termination of pregnancy.
- Shaped distal end with two large lateral eyes.
- Size : 4.0 to12mm.

M-7610
SURGICAL GOWN



- Gown material available Non-Woven & SMS.
- GSM : 43GSM.
- Sterile / Non - Sterile Both Available.

M-7607
DISPOSABLE CAP

BOUFFANT CAP



SURGEON CAP



M-7604
SURGICAL GLOVES (STERILE)



- Made from pure natural rubber latex, powdered with biologically absorbable powder.
- Extremely soft rubber reduces finger fatigue and increases efficiency.

M-7606
SURGICAL FACE MASK(3PLY)



- Mask with melt blown filter.
- Flexible ear loop.

M-7611
N95 FACE MASK



- 5-PLY Mask with three / two melt blown filter.
- Adjustable nose clip and flexible ear loop.

M-7616
DISPOSABLE SHOE COVER



- Ankle length.
- Non-woven.

Authorised Marketing & Distribution
(GCC, MENA, AFRICA)



P.O.Box: 62283, WH-2, 10C Street,
Ras Al Khor Industrial Area 1, Dubai, UAE.
E-mail : contact@justcaremed.com

www.umaflowindia.com



Plot no. 832, Nr. Ganesh Rubber,
Prima Automization Lane, Santej, Tal. Kalol,
Dist. Gandhinagar, Gujarat, INDIA.
M : +91 97243 12823
E : royalsurgicare@gmail.com
International Inquiry  +91 97268 65503