

Hospital Furniture 2019-20

KRAFT®
Crafted for Care & Comfort



 **Justcare**

ICU Beds



KRAFT 100
MOTORIZED ICU BED
5 FUNCTION (EXCEL PLUS)



KRAFT 101
MOTORIZED ICU BED
5 FUNCTION (EXCEL)



KRAFT 102
MOTORIZED ICU BED
5 FUNCTION (PREMIUM)



KRAFT 103
MANUAL ICU BED
5 FUNCTION (DELUXE)



KRAFT 121
MOTORIZED ICU BED
3 FUNCTION (PREMIUM)



KRAFT 122
MANUAL ICU BED
3 FUNCTION (DELUXE)

Product Code	Type	Functions	Bed surface	Head & Leg Bows	Side Rails	4 No. of Wheels	Powder Coating	Includes/Optional
KRAFT 100	Motorized	5	P.P. Bed	P.P	4 P.P Tuck Away	125mm - all with central locking system	9 Tank Anti microbial	Includes: CPR, battery back up and nursing panel, IV Rod
KRAFT 101	Motorized	5	P.P. Bed	P.P	4 P.P Tuck Away	125mm - Twin wheels with brakes	9 Tank Anti microbial	Optional: battery backup and Nursing panel
KRAFT 102	Motorized	5	Steel Pressed	P.P	2 Steel Collaps - 6 Pole	125mm - all with brakes	9 Tank Epoxy	Optional: ABS Side Rails
KRAFT 103	Manual	5	Steel Pressed	P.P	2 Steel Collaps - 6 Pole	125mm - all with brakes	9 Tank Epoxy	Optional: ABS Side Rails
KRAFT 121	Motorized	3	Steel Pressed	P.P	2 Steel Collaps - 6 Pole	125mm - all with brakes	9 Tank Epoxy	Optional: ABS Side Rails
KRAFT 122	Manual	3	Steel Pressed	P.P	2 Steel Collaps - 6 Pole	125mm - all with brakes	9 Tank Epoxy	Optional: ABS Side Rails

Abbreviation: P.P - Polypropylene & M.S - Mild Steel

5 FUNCTION



3 FUNCTION



Ward Care Beds



KRAFT 131
MOTORIZED FOWLER BED
2 FUNCTION (PREMIUM)



KRAFT 132
MANUAL FOWLER BED
2 FUNCTION (DELUXE)



KRAFT 133
MANUAL FOWLER BED
2 FUNCTION (ECONOMY)



KRAFT 134
MANUAL FOWLER BED
2 FUNCTION (CLASSIC)



KRAFT 141
MOTORIZED SEMI FOWLER BACKREST BED
1 FUNCTION (PREMIUM)



KRAFT 142
MANUAL SEMI FOWLER BACKREST BED
1 FUNCTION (DELUXE)



KRAFT 143
MANUAL SEMI FOWLER BACKREST BED
1 FUNCTION (ECONOMY)



KRAFT 144
MANUAL SEMI FOWLER BACKREST BED
1 FUNCTION (CLASSIC)



KRAFT 123
MANUAL ICU BED
3 FUNCTION (ECONOMY)

Product Code	Type	Functions	Bed surface	Head & Leg Bows	Side Rails (optional)	4 No. of Wheels(optional)
KRAFT 131	Motorized	2	Steel Pressed	P.P	2 Steel Collaps - 6 Pole	125mm - all with brakes
KRAFT 132	Manual	2	Steel Pressed	P.P	2 Steel Collaps - 6 Pole	125mm - all with brakes
KRAFT 133	Manual	2	Strip model	P.P	2 Alumi Collaps - 5 Pole	125mm - all with brakes
KRAFT 134	Manual	2	Strip model	M.S	2 M.S. Collapsible	Nil
KRAFT 141	Motorized	1	Steel Pressed	P.P	2 Steel Collaps - 6 Pole	125mm - all with brakes
KRAFT 142	Manual	1	Steel Pressed	P.P	2 Steel Collaps - 6 Pole	125mm - all with brakes
KRAFT 143	Manual	1	Strip model	P.P	2 Alumi Collaps - 5 Pole	125mm - all with brakes
KRAFT 144	Manual	1	Strip model	M.S	2 M.S. Collapsible	Nil
KRAFT 123	Manual	3	Strip model	P.P	2 Alumi Collaps - 5 Pole	125mm - all with brakes

2 FUNCTION



Back Rest Lifting



Leg Rest Lifting

1 FUNCTION



Back Rest Lifting



KRAFT 151
PLAIN BED CLASSIC
SIZE : L 1900mm X W 900 mm X H 560mm



KRAFT 152
ATTENDANT BED
SIZE : L 1830mm X W 610mm X H 380mm



KRAFT 153
PEDIATRIC BED CLASSIC
SIZE : L 1370mm X W 625mm X H 560mm



KRAFT 154
ATTENDANT BED CUM CHAIR



KRAFT 201
EXAMINATION COUCH
SIZE : L 1830mm X W 575 mm X H 850mm



KRAFT 202
EXAMINATION COUCH WITH BACKREST
SIZE : L 1830mm X W 575 mm X H 850mm



KRAFT 203
SIMPLE EXAMINATION COUCH
SIZE : L 1830mm X W 575mm X H 850mm



KRAFT 204
EXAMINATION TABLE S.S.
SIZE : L 1830mm X W 537mm X H 810mm



KRAFT 205
BLOOD DONATION COUCH



KRAFT 301
MOTORIZED OBSTETRIC LABOR &
RECOVERY BED - 3 FUNCTION



KRAFT 302
TELESCOPIC OBSTETRIC & GYNAEC TABLE - S.S.
SIZE : L 1830mm X W 760 mm X H 760mm



KRAFT 303
SIMPLE GYNAEC DELIVERY TABLE S.S./M.S.
SIZE : L 1830mm X W 760mm X H 760mm



KRAFT 400
TRAUMA CARE RECOVERY TROLLEY
By Screw Mecahnism

KRAFT 401
HYDRAULIC TRAUMA CARE
RECOVERY TROLLEY

KRAFT 402
STRETCHER ON TROLLEY - SS
SIZE : L 2030mm X W 610mm X H 870mm

KRAFT 403
ABS PATIENT TRANSFER TROLLEY

KRAFT 404
ALUMINIUM FOLDING STRETCHER



KRAFT 500 A
ABS EMERGENCY TROLLEY

KRAFT 500 B
ABS ANESTHESIA TROLLEY

KRAFT 501
EMERGENCY CRASH CART - S.S

KRAFT 502
BABY TROLLEY - S.S
SIZE : L 750mm X W 455mm X H 1000mm

KRAFT 503
MONITOR TROLLEY M.S

KRAFT 504
WASTE CARRYING TROLLEY S.S



KRAFT 505
BOWL STAND SINGLE - S.S
SIZE : H 840 mm X D 355 mm

KRAFT 506
BOWL STAND DOUBLE - S.S
SIZE : H 840 mm X D 355 mm

KRAFT 507
REVOLVING STOOL - S.S
SIZE : W 290mm X H 425mm to 675mm

KRAFT 508
SOILED LINEN TROLLEY BIG - S.S
SIZE : L 790mm W 560mm X H 1030mm

KRAFT 509
BED SIDE STOOL WITH S.S TOP
SIZE : L 405mm X W 405 mm X H 455mm

KRAFT 510 / 511
SINGLE STEP STOOL - S.S / M.S
SIZE : L 450mm X B 225mm X H 225mm
DOUBLE STEP STOOL - S.S / M.S
SIZE : L 450mm X B 400mm X H 400mm

Ward / OT Furniture



KRAFT 512
SALINE STAND - S.S



KRAFT 601
DRESSING TROLLEY
Available in various sizes



KRAFT 602
INSTRUMENT TROLLEY - S.S
Available in various sizes



KRAFT 603
BED SIDE TROLLEY - S.S
Available in various sizes



KRAFT 604
INSTRUMENT TROLLEY - S.S
Available in various sizes

S.S. TROLLEY - VARIOUS SIZES

L	X	W	X	H
16	X	16	X	35
18	X	18	X	35
24	X	18	X	35
27	X	18	X	35
30	X	20	X	35
36	X	20	X	35
36	X	24	X	35
42	X	24	X	35



KRAFT 605
MAYO'S TROLLEY - S.S
SIZE : L 560mm X W 445mm X H 810 to 1270mm



KRAFT 606
OVER BED TABLE - S.S
SIZE : L 1115mm X W 445mm X H 965mm



KRAFT 701
BED SIDE LOCKER ABS (DELUXE)
SIZE : L 460mm X W 450mm X H 820mm



KRAFT 702
BED SIDE LOCKER ABS
SIZE : L 445mm X W 415mm X H 770mm



KRAFT 703 A ○ ●
ADJUSTABLE BED SIDE TABLE
GAS SPRING



KRAFT 703 B
ADJUSTABLE BED SIDE TABLE
WITH KNOB



KRAFT 704
NON FOLDING WHEEL CHAIR S.S



KRAFT 704 A
FOLDING WHEEL CHAIR



KRAFT 705
BED SIDE SCREEN S.S / M.S
SIZE : H 1675mm X W 2440mm



KRAFT 706
CYLINDER TROLLEY
S.S / M.S 10 LTRS

Accessories



KRAFT 801 ●●
P.P HEAD & LEG BOWS (ECONOMY)
approx wt: 4.500 kgs per set



KRAFT 802 ●●
P.P HEAD & LEG BOWS (DELUXE)
approx wt: 6.500 kgs per set



KRAFT 803 ●
P.P HEAD & LEG BOWS (DELUXE PLUS)
approx wt: 6.500 kgs per set



KRAFT 804 ●●●
P.P HEAD & LEG BOWS (PREMIUM)
approx wt: 7.100 kgs per set




KRAFT 805 ●
P.P HEAD & LEG BOWS (EXCEL)
approx wt: 9.000 kgs per set




KRAFT 806 ●
P.P HEAD & LEG BOWS (EXCEL PLUS)
approx wt: 9.200 kgs per set



KRAFT 807 ○●●●
TUCK AWAY P.P SIDE RAILS
(SET OF 4)



KRAFT 808 ●
TUCK AWAY P.P SIDE RAILS
(SET OF 2)



KRAFT 809 (6 Pole) / **809 A** (5 Pole)
STEEL COLLAPSIBLE SIDE RAILS
(SET OF 2)



KRAFT 810 (6 Pole) / **810 A** (5 Pole)
ALUMINIUM COLLAPSIBLE SIDE RAILS
(SET OF 2)



KRAFT 811 (3") / **811 A** (4")
SINGLE FOLD FOAM MATTRESS



KRAFT 812 (3") / **812 A** (4")
MULTI FOLD FOAM MATTRESS



KRAFT 813
CRASH CART ACCESSORIES
2 CABINET & 6 UTILITY BINS



KRAFT 814
BED CRANK - 48"



KRAFT 815
ABS EXPANDABLE FOOD TRAY



KRAFT 816
CENTRAL LOCKING CASTORS
125MM



KRAFT 817
TWIN CASTORS
125MM



KRAFT 818
HEAVY CASTORS
125MM



KRAFT 819
ECONOMY CASTORS
125MM



KRAFT 820
ABS CHART HOLDER



KRAFT 821
ABS BED SIDE TABLE TOP



KRAFT 822
I.V ATTACHMENT
FOR BED



KRAFT®
Crafted for Care & Comfort



Justcare Medical Supplies LLC

📍 P.O.Box: 62283 WH-2, 10C Street, Ras Al Khor Industrial Area 1, Dubai, UAE.
Telephone: +971-4-2869111 | Fax: +971-4-2869222 | Mobile: +971-56-4338644
Email: contact@rainstone.co, justcare.uae@gmail.com

Note:

The Products displayed in this catalog are subjected to change in design, colour, size and features due to continuous development. Customization of design, colour, size and features is offered on request. Please obtain OFFER / QUOTE with detailed specification.